



Decontamination Management of flexible endoscope decontamination (to NHS guidance)

Our **decontamination** training providers learners with hands-on experience across a wide range of equipment within our specialist **decontamination** facility.

Is this course for you?

This course provides the knowledge and skills for people involved in the management of flexible endoscope decontamination, either in dedicated endoscope units or sterile services departments; Users, decontamination leads and estate managers.

This course does NOT cover the hands-on practice of decontamination (by the Operator) and periodic testing of automated endoscope reprocessors.

Course aim

To provide learners with the necessary knowledge and skills to successfully decontaminate flexible endoscopes in a safe and effective way, and in line with current guidelines.

Objectives

- Describe the construction, care, handling and use of flexible endoscopes used in medical practice.
- Explain the process of flexible endoscope decontamination and storage.
- Explain the role of washer disinfectors in endoscope reprocessing.
- Interpret machine data for product release.
- Describe the periodic testing of flexible endoscope washer disinfectors and interpret test reports.
- Describe traceability systems for endoscope reprocessing.

Pre-requisites

Delegates attending the course should be working in an endoscopy environment and have some operational knowledge of endoscope decontamination.

50+ years of world-class training