



Decontamination SSD supervisors / managers (to NHS guidance)

Our **decontamination** training providers learners with hands-on experience across a wide range of equipment within our specialist **decontamination** facility.

Is this course for you?

This course is aimed at National Health Service staff and others in the private sector who have responsibility for the management and production of all sterile services and goods. These include Sterile Services managers, supervisors and those designated the 'User' as defined in NHS guidance.

Course aim

To enable learners to meet their responsibilities as users.

Objectives

- Define various sterilization and washer-disinfection processes
- Describe the operating principles of porous load sterilizers
- Carry out daily and weekly tests
- Identify appropriate documentation
- Explain roles and responsibilities of key personnel involved with sterile services
- Reference and describe specified testing procedures
- Correctly interpret test results provided by the CP(D) and explain certification responsibilities

Pre-requisites

Delegates should be working in or towards a supervisor/manager role within an acute trust/sterile services environment.

Additional information

The course includes practical work on sterilizers and washer-disinfectors to help in the attainment of the objectives.

Related courses

Although our Users' course covers washer disinfectors, it does not cover the more specific requirements of flexible endoscope washer disinfectors. Knowledge of this subject can be gained by attending our Management of flexible endoscope decontamination (to NHS guidance) course.

50+ years of world-class training